


# Arcadius<sup>®</sup>XP C Spinal System

Ordering Guide



Aesculap Spine

# Arcadius<sup>®XP</sup> C Spinal System

## Plasmapore<sup>®XP</sup> Enhanced ACDF Stand-Alone Interbody

### Indications for Use

The Arcadius<sup>XP</sup> C Spinal System is intended to be used as an intervertebral body fusion device as a stand-alone system used with the supplied bone screws and requires no additional supplementary fixation. It is intended for spinal fusion procedures at one level in the cervical spine from the C2-C3 disc space to the C7-T1 disc space for the treatment of cervical degenerative disc disease (defined as discogenic pain with degeneration of the disc confirmed by history and radiographic studies) using autograft bone. Patients should be skeletally mature and must have undergone a regimen of at least six (6) weeks of nonoperative treatment prior to being treated with the Arcadius<sup>XP</sup> C Spinal System.

### Contraindications

Any medical or surgical condition that could preclude the potential success of the implantation. These include:

- Existing or risk of acute or chronic infections, fever/leukocytosis
- Progressive joint disease or bone absorption syndromes such as Paget's disease, osteopenia, osteoporosis or osteomyelitis which may prevent adequate fixation. Severe defects of the bony structures of the spine which are a prerequisite for stable implantation of the cages.
- Bone tumors in the region of the implant anchoring
- Proven or suspected allergy/foreign body reactions to implant materials
- Conditions that may place excessive stresses on bone and implants, such as severe obesity, pregnancy or degenerative diseases. The decision to use this system in such conditions must be made by the physician, taking into account the risks versus the benefits to the patient.
- Patients resistant to following postoperative restrictions on movement, especially in athletic and occupational activities
- Systemic, metabolic and degenerative diseases
- Psychosocial problems; unwillingness or inability of the patient to follow the instructions for postoperative treatment
- Drug abuse or alcoholism
- Generally poor condition of the patient
- Prior fusion at the level(s) to be treated
- All cases that are not listed under "Indications"

For a complete list of Warnings, Precautions and Risks, please visit our website at [www.aesculapimplantsystems.com](http://www.aesculapimplantsystems.com).

### Features


- Wide variety of implant options
- Generous graft window
- Surface texturing
- Two X-ray marker pins



Self-tapping, self-drilling 4 mm diameter bone screws

Diverging screw design with 35 degree cranial-caudal orientation

- Diverging screw design
- Dual locking mechanism
- Self-tapping bone screws
- Comprehensive array of instrumentation



Plasmapore<sup>XP</sup>  
Surface Technology



Screw integrates with PEEK  
threading and internal locking rim  
to create dual locking mechanism

## Ordering Information

Implants are delivered sterile packed.

### Set Numbers

**ST0632** Arcadius<sup>XP</sup> C Implant Bank

**ST0652** Arcadius<sup>XP</sup> C Instrument Set

### Implant Bank Configuration

The Arcadius<sup>XP</sup> C Implant Set has been optimized to include the following quantities of each implant footprint (Height x Length in mm):

#### Implant Bank Configuration (ST0632) 13 x 16 mm Footprint Width

Height	Quantity	4 Degrees	7 Degrees
5 mm	3	S0706P	S0726P
6 mm	3	S0707P	S0727P
7 mm	3	S0708P	S0728P
8 mm	2	S0709P	S0729P
9 mm	2	S0710P	S0730P
10 mm	1	S0711P	S0731P
11 mm	1	S0712P	S0732P

#### Implant Bank Configuration (ST0632) 15 x 17 mm Footprint Width

Height	Quantity	4 Degrees	7 Degrees
5 mm	3	S0746P	S0766P
6 mm	3	S0747P	S0767P
7 mm	3	S0748P	S0768P
8 mm	2	S0749P	S0769P
9 mm	2	S0750P	S0770P
10 mm	1	S0751P	S0771P
11 mm	1	S0752P	S0772P

### Bone Screws

Item No.	Quantity	Description	Color
S0791T	8	14 mm Bone Screw	Blue
S0792T	8	16 mm Bone Screw	Gold
S0793T	8	18 mm Bone Screw	Green

- Stand-alone, **zero-profile** design avoids supplementary fixation
- Dual **locking mechanism** to help prevent screw back-out
- Anatomical implant with **28 sizes** for individualized patient care
- Variations of **key instruments** to accommodate surgeon preference
- Diverging screws **prevent interference** with cervical plate screws in adjacent level constructs

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2018 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA.  
Aesculap is an equal opportunity employer

Aesculap Implant Systems, LLC | 3773 Corporate Parkway | Center Valley, PA | 18034  
Phone 866-229-3002 | Fax 610-984-9096 | [www.aesculapimplantsystems.com](http://www.aesculapimplantsystems.com)

Aesculap Implant Systems, LLC - a B. Braun company

DOC1511 Rev A 500 2/18